

MINI PLUS CART



- Designed to bring dressing set, instruments and medicine with easy when doctors and nurses for treating the patients
- All stainless steel 304 construction with non spills edge and handle
- Well stabilized and strong structured (use corner ring)
- Calmed and good rolling 4" casters with brake

● SPECIFICATION

Body	1.2T SUS 304/27 plate
Shelf	1.2T SUS 304/27 plate
Drawer	SUS 304
Leg	Ø25.4x1.5T SUS 304/27 pipe
Guard rail	Ø9.6 SUS 304/27 bar
Cornering	Die casting
Bumper	Ø90 round rubber bumper
Caster	4 inch

● Option Accessories



- Handle
Ergonomic designed handle
- Guard rail
guard rail for not dropping of fluid or instrument



- Shelf
About 15mm indented top tray helps to prevent things falling off



- Disposal can
Size : φ 120x140mm
φ 200x210mm
φ 250x260mm



- Caster
Calmed and good rolling casters with brake



- Shelf for waste disposal box